





UX Version: ADI_UX_210105_V2.1

Designer :Holiday




Starting up、Restart、Uparade

starting up




In power off mode, press the Navigation button for 1s to start the machine

starting up in shipment mode




It must be started with a charging cable


restart




upgrade




Main interface (static page)





 short press: Navigation function


No action in current Page: Backlight will auto close after 30s

Heart rate




-  short press: Confirm
-  long press: Quit

Waveform




No action in current Page: Backlight will auto close after 30s

LOW_TOUCH LOGGING




 short press: Confirm




 long press: Quit/Back

 short press: Navigation function

No action in current Page: Backlight will auto close after 30s

Setting (static page)



-  short press: Confirm
-  long press: Quit/Back
-  short press: Navigation function

No action in current Page: Backlight will auto close after 30s