

Machine to Machine Messaging V.2 (M2M2) Commands

RAJESH VIJAYAKUMAR



1/18/2021

M2M2 Header

- Based on variable-length UDP packets (2 bytes each for source, destination, checksum, length)
- Application-focused
 - All applications have an address
 - Packets are dynamically routed between applications and across machines based on addresses
- Implements Point-to-Point, Publish/Subscribe, and Broadcast messaging
 - A "Publisher" of data doesn't have to know who is consuming that data, or where they are
- ► All inter-task messaging is done over M2M2
 - Individual tasks can be addressed transparently from anywhere in the system
 - Allows applications to be easily moved between machines in the system

0 1	2 3	4 5	6 7	8 9
Source	Destination	Length	Checksum	Payload



M2M2 – Commands Classification

- Each application or stream has an address This is used to route packets across applications or machines
- Various sensor/bio-medical applications/streams supported
 - ADPD ADPD stream
 - ADXL ADXL stream
 - PPG PPG stream;
 - ECG ECG stream
 - EDA EDA stream
 - Temperature temperature stream
 - Pedometer pedometer stream
 - Synchronized ADPD and ADXL sync ADPD/ADXL stream
 - AD5940
- House-keeping/Utility applications/streams
 - File System
 - System application Battery Stream information




M2M2 – Getting Started with the sensor/bio-medical application

- Define the application and its stream address
 - M2M2_ADDR_MED_ECG = 49921 → Application
 - M2M2_ADDR_MED_ECG_STREAM = 50177 → Stream
- Define an interface file for the application Interface contains all structs and enums used to interact with an application
 - Eg:- ecg_application_interface.py single source of truth
 - The common application specific commands can be inherited from common_application_interface
 - The common sensor specific commands can be inherited from common_sensor_interface
 - Using data_definition_generator convert into language-specific definitions C, C++ and python
 - Provides a single place to update the interfaces and have them propagate through the entire system (and tools)
- Application Main components will be an Init(), message send() and process().
 - M2M2_APP_COMMON_CMD_STREAM_START_REQ → Init() → M2M2_APP_COMMON_CMD_STREAM_START_RESP
 - M2M2_APP_COMMON_CMD_STREAM_SUBSCRIBE_REQ → M2M2_APP_COMMON_CMD_STREAM_SUBSCRIBE_RESP
 - Process() process data received from mailbox
 - Message_send() → All command responses
 - M2M2_APP_COMMON_CMD_STREAM_STOP_REQ → Deinit() → M2M2_APP_COMMON_CMD_STREAM_STOP_RESP
 - M2M2_APP_COMMON_CMD_STREAM_UNSUBSCRIBE_REQ → M2M2_APP_COMMON_CMD_STREAM_UNSUBSCRIBE_RESP




M2M2 – Getting Started with the sensor/bio-medical application.....





M2M2 – Command/Response Packets





M2M2 Tooling - Python Command Line Interface (CLI)

- Written in Python
- Just an M2M2 application
- Used as the reference tool for development
 - Firmware applications are developed and tested with the CLI first
- Easy to add new commands
- Easy to extend
- Easy to automate
- Built-in helptext

```
OTE: Not connected to a serial device!
et status of a sensor or application.
vailable devices:
    : The ECG service.
     The EDA service.
      The PPG heart rate service.
      The ADXL device.
 yncppg': The sync PPG data stream service.
 dpd': The ADPD device.
   #>status [device]
   #>status adpd
#>status adxl
   #>status ppg
OTE: Not connected to a serial device!
Set/Get the status of the ADP5350's DCDC converter.
his has the effect of enabling or disabling USB bus power to the device.
   #>usbPwr [get/enable/disable]
   #>usbPwr enable
```

```
C:\Users\jzahn\Documents\git\gen3\m2M2\tools>CLI.py
This is the m2m2 UART shell. Type "help" or "?" to list commands.
NOTE: Not connected to a serial device!
Documented commands (type help <topic>):
batteryTest
                    fs_vol_info
                                                  getVersion
clockCálibration getAdpdVersion
                                                                  sensor
connect
                    getAdxlVersion
                                                                 setBatteryCharging
setBatteryThreshold
                    getBatteryInfo
flush
                    getDateTime
fs_format
                    getDcfg
                                                                 setDateTimePS
fs_log
fs_ls
                    getEdaVersion
                                                                  setLed
                    aetLcfa
                                                                 setPowerMode
                                                  loadDevice
                                                                 setPowerModePS
fs_mount
                    getLed
 s_refhr
                    getPpgAlgoVendorVersion
                                                 msg_verbose
                                                                 setPpgLcfg
                    getPpgStateInfo
                                                                  setPpgSync
setSlot
 s_rm
 s_status
                    getPpgStates
                                                  quickstart
 s_stream
                    getPpgSync
                                                                 status
                                                  quickstop
 s_sub
 s_sub_status
                                                                  usbPwr
Undocumented commands:
```



M2M2 Tooling - Basic CLI Commands for an application

- getVersion
 - Gets Version of PM and PS
- fs_log_start/stop
 - Start or stop logging
- fs_stream <file>
 - Reads a file from the NAND flash
 - Tracks file transfer rate and progress
- quickstart ecg
 - Run a pre-defined set of CLI commands with one shortcut
- msg_verbose
 - Changes the CLI's verbosity level
 - Can dump raw hex values of packets
 - Coloured messages to distinguish between verbosity levels

