



UX Version: ADI_UX_210304_V2.9

Designer :Holiday



Action button (Setting)



Navigation button (Function)

Gestures specifications

Function page	Navigation button		Navigation function
	Action button		Confirm
	Action button Long press		Quit/Back

In shutdown mode	Navigation button		Starting up
work mode	Navigation button Long press		Restart
work mode	Navigation button + Action button Long press		Upgrade
Shipment mode	USB		Charging activation
			Automatic jump

Dynamic Page


updating

Heart rate measurement

Waveform measurement


Starting up、Restart、Uparade

starting up




In power off mode, press the Navigation button for 1s to start the machine

starting up in shipment mode




It must be started with a charging cable


restart



upgrade




Main interface (static page)



short press: Navigation function

No action in current Page: Backlight will auto close after 30s

Heart rate




short press: Confirm




long press: Quit




Waveform



No action in current Page: Backlight will auto close after 30s


LOW_TOUCH LOGGING






-  short press: Confirm
-  long press: Quit/Back
-  short press: Navigation function

No action in current Page: Backlight will auto close after 30s


Setting (static page)



-  short press: Confirm
-  long press: Quit/Back
-  short press: Navigation function


No action in current Page: Backlight will auto close after 30s

ABOUT



 short press: Confirm

 long press: Quit/Back

 short press: Navigation function

No action in current Page: Backlight will auto close after 30s